[suPAR]

Monografía

Contenido:

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1. Ploug M1, Eriksen J, Plesner T, Hansen NE, Danø K.A soluble form of the glycolipid-anchored receptor for urokinase-type plasminogen activator is secreted from peripheral blood leukocytes from patients with paroxysmal nocturnal hemoglobinuria. *Eur J Biochem.* 1992 Sep 1;208(2):397-404.
2. Eugen-Olsen, J. et al. Circulating soluble urokinase plasminogen activator receptor predicts cancer, cardiovascular disease, diabetes and mortality in the general population. *J. Intern. Med.* 2010;268, 296–308.
3. Borne Y, Persson M, Melander O, et al. Increased plasma level of soluble urokinase plasminogen activator receptor is associated with incidence of heart failure but not atrial fibrillation. *European journal of heart failure* 2014;16:377-83.
4. Eapen DJ, Manocha P, Ghasemzedah N, et al. Soluble urokinase plasminogen activator receptor level is an independent predictor of the presence and severity of coronary artery disease and of future adverse events. *Journal of the American Heart Association* 2014;3:e001118.
5. Gumus A, Altintas N, Cinarka H, et al. Soluble urokinase-type plasminogen activator receptor is a novel biomarker predicting acute exacerbation in COPD. *International journal of chronic obstructive pulmonary disease* 2015;10:357-65.

6. Lyngbaek S, Andersson C, Marott JL, et al. Soluble urokinase plasminogen activator receptor for risk prediction in patients admitted with acute chest pain. *Clinical chemistry* 2013;59:1621-9.
7. Lyngbaek S, Marott JL, Moller DV, et al. Usefulness of soluble urokinase plasminogen activator receptor to predict repeat myocardial infarction and mortality in patients with ST-segment elevation myocardial infarction undergoing primary percutaneous intervention. *The American journal of cardiology* 2012;110:1756-63.
8. Molkanen T, Ruotsalainen E, Thorball CW, et al. Elevated soluble urokinase plasminogen activator receptor (suPAR) predicts mortality in *Staphylococcus aureus* bacteremia. *Eur J Clin Microbiol Infect Dis*. 2011;30:1417-24.
9. Persson M, Ostling G, Smith G, et al. Soluble urokinase plasminogen activator receptor: a risk factor for carotid plaque, stroke, and coronary artery disease. *Stroke; a journal of cerebral circulation* 2014;45:18-23.
10. Raggam RB, Wagner J, Pruller F, et al. Soluble urokinase plasminogen activator receptor predicts mortality in patients with systemic inflammatory response syndrome. *Journal of internal medicine* 2014;276:651-8.
11. Hayek S, Sever S, Ko Y et al. Soluble urokinase receptor and chronic kidney disease. *N Engl J Med*. 2015 Nov 12; 373(20):1916-25.
12. Reichsoellner M, Raggam RB, Wagner J, Krause R, Hoenigl M. Clinical evaluation of multiple inflammation biomarkers for diagnosis and prognosis for patients with systemic inflammatory response syndrome. *Journal of clinical microbiology* 2014;52:4063-6.
13. Suberviola B, Castellanos-Ortega A, Ruiz Ruiz A, Lopez-Hoyos M, Santibanez M. Hospital mortality prognostication in sepsis using the new biomarkers suPAR and proADM in a single determination on ICU admission. *Intensive care medicine* 2013;39:1945-52.
14. Zimmermann HW, Koch A, Seidler S, Trautwein C, Tacke F. Circulating soluble urokinase plasminogen activator is elevated in patients with chronic liver disease, discriminates stage and aetiology of cirrhosis and predicts prognosis. *Liver international : official journal of the International Association for the Study of the Liver* 2012;32:500-9.
15. Rasmussen LJH, Ladelund S, Haupt TH, Ellekilde GE, Eugen-Olsen J, Andersen O. Combining National Early Warning Score With Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. *Crit Care Med*. 2018 Dec;46(12):1961-1968.
16. Donadello K, Scolletta S, Taccone FS, et al. Soluble urokinase-type plasminogen activator receptor as a prognostic biomarker in critically ill patients. *Journal of critical care* 2014;29:144-9.
17. Lyngbaek S, Marott JL, Sehestedt T, et al. Cardiovascular risk prediction in the general population with use of suPAR, CRP, and Framingham Risk Score. *International journal of cardiology* 2013;167:2904-11.
18. Sidenius N, Sier CF, Ullum H et al. Serum level of soluble urokinase-type plasminogen activator receptor is a strong and independent predictor of survival in human immunodeficiency virus infection. *Blood* 2000;96(13):4091-4095.
19. Eugen-Olsen J, Gustafson P, Sidenius N et al. The serum level of soluble urokinase receptor is elevated in tuberculosis patients and predicts mortality during treatment: a community study from Guinea-Bissau. *Int J Tuberc Lung Dis* 2002;6(8):686-692.
20. Perch M, Kofoed P, Fischer TK et al. Serum levels of soluble urokinase plasminogen activator receptor is associated with parasitemia in children with acute *Plasmodium falciparum* malaria infection. *Parasite Immunol* 2004;26(5):207-211.
21. Tzanakaki G, Pappas M, Kyprianou M, Barbouni A, Eugen-Olsen J, Kourea-Kremastinou J. Elevated soluble urokinase receptor values in CSF, age and bacterial meningitis infection are independent and additive risk factors of fatal outcome. *Eur J Clin Microbiol Infect Dis* 2012;31(6):1157-1162.
22. Wrotek A, Jackowska T. The role of the soluble urokinase plasminogen activator (suPAR) in children with pneumonia. *Respir Physiol Neurobiol* 2015;209:120-3. doi: 10.1016/j.resp.2014.12.018. Epub: 2015 Jan 17.:120-123.
23. Andersen O, Eugen-Olsen J, Kofoed K, Iversen J, Haugaard SB. Soluble urokinase plasminogen activator receptor is a marker of dysmetabolism in HIV-infected patients receiving highly active antiretroviral therapy. *J Med Virol*. 2008 Feb;80(2):209-16.

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-
1. Eugen-Olsen, J. *et al.* Circulating soluble urokinase plasminogen activator receptor predicts cancer, cardiovascular disease, diabetes and mortality in the general population. *J. Intern. Med.* 268, 296–308 (2010).
 2. Borné, Y., Persson, M., Melander, O., Smith, J. G. & Engström, G. Increased plasma level of soluble urokinase plasminogen activator receptor is associated with incidence of heart failure but not atrial fibrillation. *Eur. J. Heart Fail.* 16, 377–383 (2014).
 3. Haupt, T. H. *et al.* Risk factors associated with serum levels of the inflammatory biomarker soluble urokinase plasminogen activator receptor in a general population. *Biomed. Insights* 9, 91–100 (2014).
 4. Persson, M., Engström, G., Björkbacka, H. & Hedblad, B. Soluble urokinase plasminogen activator receptor in plasma is associated with incidence of CVD. Results from the Malmö Diet and Cancer Study. *Atherosclerosis* 220, 502–505 (2012).
 5. Langkilde, A. *et al.* Increased plasma soluble uPAR level is a risk marker of respiratory cancer in initially cancer-free individuals. *Cancer Epidemiol. Biomark. Prev. Publ. Am. Assoc. Cancer Res. Cosponsored Am. Soc. Prev. Oncol.* 20, 609–618 (2011).
 6. Persson, M., Engström, G., Björkbacka, H. & Hedblad, B. Soluble urokinase plasminogen activator receptor in plasma is associated with incidence of CVD. Results from the Malmö Diet and Cancer Study. *Atherosclerosis* 220, 502–505 (2012).
 7. Hayek, S. S. *et al.* Soluble Urokinase Receptor and Chronic Kidney Disease. *N. Engl. J. Med.* 373, 1916–1925 (2015).
 8. Haupt TH *et al.* Healthy lifestyles reduce suPAR and mortality in a Danish general population study. *Immun Ageing.* 22;16:1 (2019).
 9. Zimmermann, H. W., Koch, A., Seidler, S., Trautwein, C. & Tacke, F. Circulating soluble urokinase plasminogen activator is elevated in patients with chronic liver disease, discriminates stage and aetiology of cirrhosis and predicts prognosis. *Liver Int. Off. J. Int. Assoc. Study Liver* 32, 500–509 (2012).
 10. Törnkvist PBS *et al.* Soluble urokinase plasminogen activator receptor is linearly associated with dietary quality and predicts mortality. *Br J Nutr.* 2019;121(6):699-708.

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1. Haupt, T. H. et al. Plasma suPAR levels are associated with mortality, admission time, and Charlson Comorbidity Index in the acutely admitted medical patient: a prospective observational study. *Crit. Care* 16, R130 (2012).
2. Rasmussen, L. J. H. et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study *Emerg Med J.* 2016 Nov;33(11):769-775
3. Raggam, R. B. *et al.* Soluble urokinase plasminogen activator receptor predicts mortality in patients with systemic inflammatory response syndrome. *J. Intern. Med.* **276**, 651–8 (2014).
4. Nayak, R. K., Allingstrup, M., Phanareth, K. & Kofoed-enevoldsen, A. suPAR as a biomarker for risk of readmission and mortality in the acute medical setting. 62, 1–4 (2015).
5. Rasmussen LJH et al, Combining National Early Warning Score With Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. *Crit Care Med.* 2018 Dec;46(12):1961-1968
6. Schultz M at al. Availability of suPAR in emergency departments may improve risk stratification: a secondary analysis of the TRIAGE III trial. *Scand J Trauma Resusc Emerg Med.* 2019 Apr 11;27(1):43. Martin Schultz, Line J. H. Rasmussen, Thomas Høi-Hansen, et al., "Early Discharge from the
7. Schultz et al. Emergency Department Based on Soluble Urokinase Plasminogen Activator Receptor (suPAR) Levels: A TRIAGE III Substudy," *Disease Markers*, vol. 2019, Article ID 3403549, 8 pages, 2019

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-
1. Sidenius N, Sier CF, Ullum H et al. Serum level of soluble urokinase-type plasminogen activator receptor is a strong and independent predictor of survival in human immunodeficiency virus infection. *Blood* 2000;96(13):4091-4095.
 2. Ostrowski SR, Katzenstein TL, Pedersen M et al. Plasma levels of intact and cleaved urokinase receptor decrease in HIV-1-infected patients initiating highly active antiretroviral therapy. *Scand J Immunol* 2006;63(6):478-486.

3. Andersen O, Eugen-Olsen J, Kofoed K, Iversen J, Haugaard SB. Soluble urokinase plasminogen activator receptor is a marker of dysmetabolism in HIV-infected patients receiving highly active antiretroviral therapy. *J Med Virol* 2008;80(2):209-216.
4. Langkilde A, Petersen J, Klausen HH, Henriksen JH, Eugen-Olsen J, Andersen O. Inflammation in HIV-infected patients: impact of HIV, lifestyle, body composition, and demography - a cross sectional cohort study. *PLoS One* 2012;7(12):e51698.
5. Hoenigl M et al and the ACTG NWCS 411 study team. Soluble Urokinase Plasminogen Activator Receptor (suPAR) is predictive of Non-AIDS Events during Antiretroviral Therapy-mediated Viral Suppression. *Clin Infect Dis*. 2018 ciy966
6. Donadello K, Scolletta S, Covajes C, Vincent JL. suPAR as a prognostic biomarker in sepsis. *BMC Med* 2012;10:2. doi: 10.1186/1741-7015-10-2.:2-10.
7. Donadello K, Scolletta S, Taccone FS et al. Soluble urokinase-type plasminogen activator receptor as a prognostic biomarker in critically ill patients. *J Crit Care* 2014;29(1):144-149.
8. Koch A, Voigt S, Kruschinski C et al. Circulating soluble urokinase plasminogen activator receptor is stably elevated during the first week of treatment in the intensive care unit and predicts mortality in critically ill patients. *Crit Care* 2011;15(1):R63.
9. Suberviola B, Castellanos-Ortega A, Ruiz RA, Lopez-Hoyos M, Santibanez M. Hospital mortality prognostication in sepsis using the new biomarkers suPAR and proADM in a single determination on ICU admission. *Intensive Care Med* 2013;39(11):1945-1952.
10. Giamarellos-Bourboulis EJ, Georgitsi M. Host response biomarker in sepsis: suPAR detection. *Methods Mol Biol* 2015;1237:241-6. doi: 10.1007/978-1-4939-1776-1_18.:241-246.
11. Giamarellos-Bourboulis EJ, Norrby-Teglund A, Mylona V et al. Risk assessment in sepsis: a new prognostication rule by APACHE II score and serum soluble urokinase plasminogen activator receptor. *Crit Care* 2012;16(4):R149.
12. Sevgi DY, Bayraktar B, Gunduz A et al. Serum soluble urokinase-type plasminogen activator receptor and interferon-gamma-induced protein 10 levels correlate with significant fibrosis in chronic hepatitis B. *Wien Klin Wochenschr* 2015.
13. Berres ML, Schlosser B, Berg T, Trautwein C, Wasmuth HE. Soluble urokinase plasminogen activator receptor is associated with progressive liver fibrosis in hepatitis C infection. *J Clin Gastroenterol* 2012;46(4):334-338.
14. Sjowall C, Martinsson K, Cardell K, Ekstedt M, Kechagias S. Soluble urokinase plasminogen activator receptor levels are associated with severity of fibrosis in nonalcoholic fatty liver disease. *Transl Res* 2015;165(6):658-666.
15. Eugen-Olsen J, Gustafson P, Sidenius N et al. The serum level of soluble urokinase receptor is elevated in tuberculosis patients and predicts mortality during treatment: a community study from Guinea-Bissau. *Int J Tuberc Lung Dis* 2002;6(8):686-692.
16. Rabna P, Andersen A, Wejse C et al. High mortality risk among individuals assumed to be TB-negative can be predicted using a simple test. *Trop Med Int Health* 2009;14(9):986-994.
17. Rabna P, Andersen A, Wejse C et al. Utility of the plasma level of suPAR in monitoring risk of mortality during TB treatment. *PLoS One* 2012;7(8):e43933.
18. Perch M, Kofoed P, Fischer TK et al. Serum levels of soluble urokinase plasminogen activator receptor is associated with parasitemia in children with acute *Plasmodium falciparum* malaria infection. *Parasite Immunol* 2004;26(5):207-211.
19. Plewes K, Royakkers AA, Hanson J et al. Correlation of biomarkers for parasite burden and immune activation with acute kidney injury in severe *falciparum* malaria. *Malar J* 2014;13:91. doi: 10.1186/1475-2875-13-91.:91-13.
20. Ostrowski SR, Ullum H, Goka BQ et al. Plasma concentrations of soluble urokinase-type plasminogen activator receptor are increased in patients with malaria and are associated with a poor clinical or a fatal outcome. *J Infect Dis* 2005;191(8):1331-1341.
21. Tzanakaki G, Paparoupa M, Kyprianou M, Barbouni A, Eugen-Olsen J, Kourea-Kremastinou J. Elevated soluble urokinase receptor values in CSF, age and bacterial meningitis infection are independent and additive risk factors of fatal outcome. *Eur J Clin Microbiol Infect Dis* 2012;31(6):1157-1162.
22. Ostergaard C, Benfield T, Lundgren JD, Eugen-Olsen J. Soluble urokinase receptor is elevated in cerebrospinal fluid from patients with purulent meningitis and is associated with fatal outcome. *Scand J Infect Dis* 2004;36(1):14-19.

23. Wrotek A, Jackowska T. The role of the soluble urokinase plasminogen activator (suPAR) in children with pneumonia. *Respir Physiol Neurobiol* 2015;209:120-3. doi: 10.1016/j.resp.2014.12.018. Epub; %2015 Jan 17.:120-123.
24. Savva A, Raftogiannis M, Baziaka F et al. Soluble urokinase plasminogen activator receptor (suPAR) for assessment of disease severity in ventilator-associated pneumonia and sepsis. *J Infect* 2011;63(5):344-350.
25. Oliveira I, Andersen A, Furtado A et al. Assessment of simple risk markers for early mortality among HIV-infected patients in Guinea-Bissau: a cohort study. *BMJ Open* 2012;2(6):e001587.

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-
1. Mustjoki, S., Sidenius, N. & Sier, C. F. M. Soluble Urokinase Receptor Levels Correlate with Number of Circulating Tumor Cells in Acute Myeloid Leukemia and Decrease Rapidly during Chemotherapy Soluble Urokinase Receptor Levels Correlate with Number of Circulating Tumor Cells in Acute Myeloid Leuke. 7126–7132 (2000).
 2. Lomholt, a F., Høyer-Hansen, G., Nielsen, H. J. & Christensen, I. J. Intact and cleaved forms of the urokinase receptor enhance discrimination of cancer from non-malignant conditions in patients presenting with symptoms related to colorectal cancer. *Br. J. Cancer* 101, 992–7 (2009).
 3. Usnarska-Zubkiewicz, L., Strutyńska-Karpińska, M., Zubkiewicz-Kucharska, A., Zarębski, P. & Grabowski, K. Soluble Urokinase-Type Plasminogen Activator Receptor and Ferritin

- Concentration in Patients with Advanced Alimentary Tract Carcinoma. Relationship to Localization, Surgical Treatment and the Stage of the Disease - Preliminary Report. *Adv. Clin. Exp. Med.* 23, 959–967 (2014).
4. Tarpgaard, L. S. *et al.* Intact and cleaved plasma soluble urokinase receptor in patients with metastatic colorectal cancer treated with oxaliplatin with or without cetuximab. *Int. J. Cancer* 1, n/a–n/a (2015).
 5. Chounta, A. *et al.* Serum soluble urokinase plasminogen activator receptor as a screening test for the early diagnosis of hepatocellular carcinoma. *Liver Int.* 35, 601–607 (2015).
 6. Jing, J. *et al.* Evaluating the value of uPAR of serum and tissue on patients with cervical cancer. *J. Clin. Lab. Anal.* 26, 16–21 (2012).
 7. Cobos, E., Jumper, C. & Lox, C. Pretreatment Determination of the Serum Urokinase Plasminogen Activator and its Soluble Receptor in Advanced Small-Cell Lung Cancer or Non-Small-Cell Lung Cancer. *Clin. Appl. Thromb.* 9, 241–246 (2003).
 8. Rubio-Jurado, B. *et al.* Circulating Levels of Urokinase-Type Plasminogen Activator Receptor and D-Dimer in Patients With Hematological Malignancies. *Clin. Lymphoma Myeloma Leuk.* 15, 621–626 (2015).
 9. Henic, E., Borgfeldt, C., Christensen, I. J., Caslén, B. & Høyer-Hansen, G. Cleaved forms of the urokinase plasminogen activator receptor in plasma have diagnostic potential and predict postoperative survival in patients with ovarian cancer. *Clin. Cancer Res.* 14, 5785–93 (2008).
 10. Sorio, C. *et al.* Elevated urinary levels of urokinase-type plasminogen activator receptor (uPAR) in pancreatic ductal adenocarcinoma identify a clinically high-risk group. *BMC Cancer* 11, 448 (2011).
 11. Miyake, H. *et al.* Elevation of serum levels of urokinase-type plasminogen activator and its receptor is associated with disease progression and prognosis in patients with prostate cancer. *Prostate* 39, 123–9 (1999).
 12. Fidan, E. *et al.* Diagnostic and prognostic significance of CA IX and suPAR in gastric cancer. *Med. Oncol.* 30, 540 (2013).
 13. Langkilde, A. *et al.* Increased plasma soluble uPAR level is a risk marker of respiratory cancer in initially cancer-free individuals. *Cancer Epidemiol. Biomarkers Prev.* 20, 609–18 (2011).
 14. Garcia-Monco, J. C., Coleman, J. L. & Benach, J. L. Soluble urokinase receptor (uPAR, CD 87) is present in serum and cerebrospinal fluid in patients with neurologic diseases. *J. Neuroimmunol.* 129, 216–223 (2002).
 15. Rasmussen L J H, *et al.* Inflammatory biomarkers and cancer: CRP and suPAR as markers of incident cancer in patients with serious nonspecific symptoms and signs of cancer. *Int J Cancer.* 2017 Jul 1;141(1):191-199.
 16. Loosen SH *et al.* Serum levels of soluble urokinase plasminogen activator receptor (suPAR) predict outcome after resection of colorectal liver metastases. *Oncotarget.* 2018 Jun 5;9(43):27027-27038.

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1. Haugaard SB, Andersen O, Hansen TW et al. The immune marker soluble urokinase plasminogen activator receptor is associated with new-onset diabetes in non-smoking women and men. *Diabet Med*. 2012 Apr;29(4):479-487
2. Heraclides A, Jensen TM, SS Rasmussen et al. The pro-inflammatory biomarker soluble urokinase plasminogen activator receptor (suPAR) is associated with incident type 2 diabetes among overweight but not obese individuals with impaired glucose regulation: effect modification by smoking and body weight status. *Diabetologia*. 2013 Apr 24
3. Theilade S, Lyngbaek S, Hansen TW, et al. Soluble urokinase plasminogen activator receptor levels are elevated and associated with complications in patients with type 1 diabetes. *J Intern Med* 2015;277
4. Rotbain Curovic V et al, Soluble Urokinase Plasminogen Activator Receptor Predicts Cardiovascular Events, Kidney Function Decline, and Mortality in Patients With Type 1 Diabetes. *Diabetes Care*. 2019. pii: dc181427.

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-
1. Haastруп E et al. Soluble urokinase plasminogen activator receptor as a marker for use of antidepressants. PLoS One. 2014 ;9(10):e110555.
 2. Ventorp F et al. Increased Soluble Urokinase-Type Plasminogen Activator Receptor (suPAR) Levels in Plasma of Suicide Attempters. PLoS One. 2015 9;10(10):e0140052.
 3. Yu L Et al. Diagnostic and prognostic significance Of suPAR in traumatic brain injury. Neurol India. 2014 Sep-Oct;62(5):498-502
 4. Rasmussen LJHH et al, Cumulative childhood risk is associated with a new measure of chronic inflammation in adulthood. J Child Psychol Psychiatry. 2019 60(2):199-208

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1. Fuhrman B. The urokinase system in the pathogenesis of atherosclerosis. *Atherosclerosis*. 2012;222(1):8-14.
2. Lyngbaek S, Marott JL, Moller DV, et al. Usefulness of soluble urokinase plasminogen activator receptor to predict repeat myocardial infarction and mortality in patients with ST-segment elevation myocardial infarction undergoing primary percutaneous intervention. *Am J Cardiol*. 2012;110(12):1756-1763.
3. Edsfeldt A, Nitulescu M, Grufman H, et al. Soluble urokinase plasminogen activator receptor is associated with inflammation in the vulnerable human atherosclerotic plaque. *Stroke*. 2012;43(12):3305-3312.
4. Blasi F, Sidenius N. The urokinase receptor: focused cell surface proteolysis, cell adhesion and signaling. *FEBS Lett*. 2010;584(9):1923-1930.
5. Hayek SS, Sever S, Ko YA, et al. Soluble Urokinase Receptor and Chronic Kidney Disease. *N Engl J Med*. 2015;373(20):1916-1925.
6. Borne Y, Persson M, Melander O, Smith JG, Engstrom G. Increased plasma level of soluble urokinase plasminogen activator receptor is associated with incidence of heart failure but not atrial fibrillation. *Eur J Heart Fail*. 2014;16(4):377-383.
7. Botha S, Fourie CM, Schutte R, Eugen-Olsen J, Pretorius R, Schutte AE. Soluble urokinase plasminogen activator receptor as a prognostic marker of all-cause and cardiovascular mortality in a black population. *Int J Cardiol*. 2015;184:631-636.
8. Sorensen MH, Gerke O, Eugen-Olsen J, et al. Soluble urokinase plasminogen activator receptor is in contrast to high-sensitive C-reactive-protein associated with coronary artery calcifications in healthy middle-aged subjects. *Atherosclerosis*. 2014;237(1):60-66.
9. Eugen-Olsen J, Andersen O, Linneberg A, et al. Circulating soluble urokinase plasminogen activator receptor predicts cancer, cardiovascular disease, diabetes and mortality in the general population. *J Intern Med*. 2010;268(3):296-308.
10. Sehestedt T, Lyngbaek S, Eugen-Olsen J, et al. Soluble urokinase plasminogen activator receptor is associated with subclinical organ damage and cardiovascular events. *Atherosclerosis*. 2011;216(1):237-243.
11. Lyngbaek S, Marott JL, Sehestedt T, et al. Cardiovascular risk prediction in the general population with use of suPAR, CRP, and Framingham Risk Score. *Int J Cardiol*. 2013;167(6):2904-2911.
12. Pawlak K, Mysliwicz M, Pawlak D. The urokinase-type plasminogen activator/its soluble receptor system is independently related to carotid atherosclerosis and associated with CC-chemokines in uraemic patients. *Thromb Res*. 2008;122(3):328-335.
13. Meijers B, Poesen R, Claes K, et al. Soluble urokinase receptor is a biomarker of cardiovascular disease in chronic kidney disease. *Kidney Int*. 2015;87(1):210-216.
14. Pawlak K, Pawlak D, Mysliwicz M. Tissue factor and urokinase-type plasminogen activator system are related to the presence of cardiovascular disease in hemodialysis patients. *Thromb Res*. 2007;120(6):871-876.
15. Intzilakis T, Hartmann G, Mouridsen MR, et al. Soluble urokinase plasminogen activator receptor, C-reactive protein and triglyceride are associated with heart rate variability in non-diabetic Danes. *Eur J Clin Invest*. 2013;43(5):457-468.

16. Lyngbaek S, Andersson C, Marott JL, et al. Soluble urokinase plasminogen activator receptor for risk prediction in patients admitted with acute chest pain. *Clin Chem*. 2013;59(11):1621-1629.
17. Botha S, Fourie CM, Schutte R, Eugen-Olsen J, Schutte AE. Soluble urokinase plasminogen activator receptor and hypertension among black South Africans after 5 years. *Hypertens Res*. 2015;38(6):439-444.
18. Rundgren M, Lyngbaek S, Fisker H, Friberg H. The inflammatory marker suPAR after cardiac arrest. *Ther Hypothermia Temp Manag*. 2015;5(2):89-94.
19. Jalkanen V, Vaahersalo J, Pettila V, et al. The predictive value of soluble urokinase plasminogen activator receptor (SuPAR) regarding 90-day mortality and 12-month neurological outcome in critically ill patients after out-of-hospital cardiac arrest. Data from the prospective FINNRESUSCI study. *Resuscitation*. 2014;85(11):1562-1567.
20. Engstrom G, Zoller B, Svensson PJ, Melander O, Persson M. Soluble urokinase plasminogen activator receptor and incidence of venous thromboembolism. *Thromb Haemost*. 2015;115(3).
21. Persson M, Ostling G, Smith G, et al. Soluble urokinase plasminogen activator receptor: a risk factor for carotid plaque, stroke, and coronary artery disease. *Stroke*. 2014;45(1):18-23.
22. Mekonnen G, Corban MT, Hung OY, et al. Plasma soluble urokinase-type plasminogen activator receptor level is independently associated with coronary microvascular function in patients with non-obstructive coronary artery disease. *Atherosclerosis*. 2015;239(1):55-60.
23. Gozdzik W, Adamik B, Gozdzik A, Rachwalik M, Kustrzycki W, Kubler A. Unchanged plasma levels of the soluble urokinase plasminogen activator receptor in elective coronary artery bypass graft surgery patients and cardiopulmonary bypass use. *PLoS One*. 2014;9(6):e98923.
24. Sörensen NA et al. Predictive value of soluble urokinase-type plasminogen activator receptor for mortality in patients with suspected myocardial infarction. *Clin Res Cardiol*. 2019 doi: 10.1007/s00392-019-01475-1.

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1. Hoenigl, M., Raggam, R.B., Wagner, J., Valentin, T., Leitner, E., Seeber, K., Zollner-Schwetz, I., Krammer, W., Prüller, F., Grisold, A.J., Krause, R., 2013. Diagnostic accuracy of soluble urokinase plasminogen activator receptor (suPAR) for prediction of bacteremia in patients with systemic inflammatory response syndrome. *Clin. Biochem.* 46, 225–229.
2. Wittenhagen, P., Kronborg, G., Weis, N., Nielsen, H., Obel, N., Pedersen, S.S., Eugen-Olsen, J., 2004. The plasma level of soluble urokinase receptor is elevated in patients with Streptococcus pneumoniae bacteraemia and predicts mortality. *Clin. Microbiol. Infect.* 10, 409–415.
3. Huttunen, R., Syrjänen, J., Vuento, R., Hurme, M., Huhtala, H., Laine, J., Pessi, T., Aittoniemi, J., 2011. Plasma level of soluble urokinase-type plasminogen activator receptor as a predictor of disease severity and case fatality in patients with bacteraemia: a prospective cohort study. *J. Intern. Med.* 270, 32–40.
4. Giamarellos-Bourboulis EJ, Norrby-Teglund A, Mylona V, Savva A, Tsangaris I, Dimopoulou I, et al. Risk assessment in sepsis: a new prognostication rule by APACHE II score and serum soluble urokinase plasminogen activator receptor. *Crit Care* 2012;16:R149.
5. Savva, A., Raftogiannis, M., Baziaka, F., Routsis, C., Antonopoulou, A., Koutoukas, P., Tsaganos, T., Kotanidou, A., Apostolidou, E., Giamarellos-Bourboulis, E.J., Dimopoulos, G., 2011. Soluble urokinase plasminogen activator receptor (suPAR) for assessment of disease severity in ventilator-associated pneumonia and sepsis. *J. Infect.* 63, 344–350.
6. Wrotek, A., Pawlik, K., Jackowska, T., 2013. Soluble receptor for urokinase plasminogen activator in community-acquired pneumonia in children. *Adv. Exp. Med. Biol.* 788, 329–334.
7. Wrotek, A., Jackowska, T., Pawlik, K., 2015. Soluble urokinase plasminogen activator receptor: an indicator of pneumonia severity in children. *Adv. Exp. Med. Biol.* 835, 1–7, http://dx.doi.org/10.1007/5584_2014_40
8. Ivancsó, I., Toldi, G., Bohács, A., Eszes, N., Müller, V., Rigó Jr., J., Vásárhelyi, B., Losonczy, G., Tamási, L., 2013. Relationship of circulating soluble urokinase plasminogen activator receptor (suPAR) levels to disease control in asthma and asthmatic pregnancy. *PLoS ONE* 8, e60697.
9. Portelli M, Siedlinski M, Stewart C et al. Genome-wide protein QTL mapping identifies human plasma kallikrein as a post-translational regulator of serum uPAR levels. *FASEB J.* 2014 Feb;28(2):923-34
10. Gumus A, Altintas N, Cinarka H, Kirbas A, Haziroglu M, Karatas M, Sahin U. Soluble urokinase-type plasminogen activator receptor is a novel biomarker predicting acute exacerbation in COPD. *Int J Chron Obstruct Pulmon Dis.* 2015 Feb 13;10:357-65.
11. Langkilde A, Hansen TW, Ladelund S, Linneberg A, Andersen O, Haugaard SB, Jeppesen J, Eugen-Olsen J Increased plasma soluble uPAR level is a risk marker of respiratory cancer in initially cancer-free individuals. *Cancer Epidemiol Biomarkers Prev.* 2011 Apr;20(4):609-18.
12. Godtfredsen NS, Jørgensen DV, Marsaa K, Ulrik CS, Andersen O, Eugen-Olsen J, Rasmussen LJH. Soluble urokinase plasminogen activator receptor predicts mortality in exacerbated COPD. *Respir Res.* 2018 21;19(1):97.

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-
1. Zimmerman HW, Koch A, Seidler S et al. Circulating soluble urokinase plasminogen activator is elevated in patients with chronic liver disease, discriminates stage and etiology of cirrhosis and predicts prognosis. *Liver Int.* 2012 Mar;32(3):500-509.
 2. Wiese S, Mortensen C, Gøtzte JO. Cardiac and proinflammatory markers predict prognosis in cirrhosis. *Liver Int.* 2014 Jul;34(6):e19-30.
 3. Tuomi, Heidi, et al. "Serum soluble urokinase plasminogen activator receptor in alcoholics: Relation to liver disease severity, fibrogenesis, and alcohol use." *J Gastroenterology Hepatic* 2014 Dec;29(12): 1991-1995.

4. Sjöwall C, Martinsson K, Cardell K et al. Soluble urokinase plasminogen activator receptor levels are associated with severity of fibrosis in nonalcoholic fatty liver disease. *Trails Res.* 2015 Jun;165(6):658-666.
5. Zimmermann HW, Reuken PA, Koch A, et al. Soluble urokinase plasminogen activator receptor is compartmentally regulated in decompensated cirrhosis and indicates immune activation and short-term mortality. *J Intern Med* 2013; 274: 86–100
Koch A, Zimmerman H, Gassler N et al. Clinical relevance and cellular source of elevated soluble urokinase plasminogen activator receptor in acute liver failure. *Liver int.* 2014 Oct;34(9):1330-1339.
6. Nikkola, A., Aittoniemi, J. J., Huttunen, R., Sand, J., & Laukkanen, J. (2015). Mo1351 Systemic Levels of Soluble Urokinase-Type Plasminogen Activator Receptor (suPAR) Predict the Severity of Acute Alcohol Pancreatitis. *Gastroenterology*, 148(4), S-680.
7. Chounta, A. *et al.* Serum soluble urokinase plasminogen activator receptor as a screening test for the early diagnosis of hepatocellular carcinoma. *Liver Int.* 35, 601–607 (2015)
8. Loosen SH et al, Serum levels of soluble urokinase plasminogen activator receptor (suPAR) predict outcome after resection of colorectal liver metastases. *Oncotarget.* 2018 5;9(43):27027-27038.

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-
1. Wei c, El Hindi S, Li J et al. Circulating urokinase receptor as a cause of focal segmental glomerulosclerosis. *Nat Med* 2011;17:952-60.
 2. Yoo H, Pedigo CE, Guzman J. Sphingomyelinase-like phosphodiesterase 3b expression levels determine podocyte injury phenotypes in glomerular disease. *J AM Soc Nephrol* 2015 jan;26(1):133-47
 3. Meijers B, Poesen R, Claes K et al. Soluble urokinase receptor is a biomarker of cardiovascular disease in chronic kidney disease. *Kidney Int* 2015 Jan;87(1):210-216.
 4. Hayek S, Sever S, Ko Y et al. Soluble urokinase receptor and chronic kidney disease. *N Engl J Med.* 2015 Nov 12; 373(20):1916-25.
 5. Musetti C, Quaglia M, Cena T et al. Circulating suPAR levels are affected by glomerular filtration rate and proteinuria in primary and secondary glomerulonephritis. *J Nephrol* 2015 Jun;28(3):299-305.
 6. Griveas I, Andriopoulos C, Sitaras P et al. Soluble urokinase plasminogen receptor: A future risk marker for hemodialysis patients *Nephrol. Dial. Transplant* 2015;30(suppl 3)
 7. Rotbain Curovic V et al, Soluble Urokinase Plasminogen Activator Receptor Predicts Cardiovascular Events, Kidney Function Decline, and Mortality in Patients With Type 1 Diabetes. *Diabetes Care.* 2019. pii: dc181427.

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-
1. Raggam RB, Wagner J, Pruller F et al. Soluble urokinase plasminogen activator receptor predicts mortality in patients with systemic inflammatory response syndrome. *J Intern Med* 2014;276(6):651-658.
 2. Koch A, Tacke F. Why high suPAR is not super-diagnostic, prognostic and potential pathogenic properties of a novel biomarker in the ICU. *Crit Care* 2011;15(6):1020.
 3. Koch A, Voigt S, Kruschinski C et al. Circulating soluble urokinase plasminogen activator receptor is stably elevated during the first week of treatment in the intensive care unit and predicts mortality in critically ill patients. *Crit Care* 2011;15(1):R63.
 4. Koch A, Zimmermann HW, Gassler N et al. Clinical relevance and cellular source of elevated soluble urokinase plasminogen activator receptor (suPAR) in acute liver failure. *Liver Int* 2014;34(9):1330-1339.
 5. Meijers B, Poesen R, Claes K et al. Soluble urokinase receptor is a biomarker of cardiovascular disease in chronic kidney disease. *Kidney Int* 2015;87(1):210-216.
 6. Giamarellos-Bourboulis EJ, Norrby-Teglund A, Mylona V et al. Risk assessment in sepsis: a new prognostication rule by APACHE II score and serum soluble urokinase plasminogen activator receptor. *Crit Care* 2012;16(4):R149.
 7. Donadello K, Scolletta S, Covajes C, Vincent JL. suPAR as a prognostic biomarker in sepsis. *BMC Med* 2012;10:2. doi: 10.1186/1741-7015-10-2-10.

8. Backes Y, van der Sluijs KF, Mackie DP et al. Usefulness of suPAR as a biological marker in patients with systemic inflammation or infection: a systematic review. *Intensive Care Med* 2012;38(9):1418-1428.
9. Eugen-Olsen J, Giamarellos-Bourboulis EJ. suPAR: The unspecific marker for disease presence, severity and prognosis. *Int J Antimicrob Agents* 2015;(15):10.
10. Backes Y, van der Sluijs KF, Tuip de Boer AM et al. Soluble urokinase-type plasminogen activator receptor levels in patients with burn injuries and inhalation trauma requiring mechanical ventilation: an observational cohort study. *Crit Care* 2011;15(6):R270.
11. Yu L, Wu X, Wang H, et al.. Diagnostic and prognostic significance of suPAR in traumatic brain injury. *Neurol India* 2014;62(5):498-502.

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-
1. Schuetz P, Hausfater P, Amin A et al. Biomarkers from distinct biological pathways improve early risk stratification in medical emergency patients: the multinational, prospective, observational TRIAGE study. *Critical Care* (2015) 19:377
 2. Biccard BM, Devereaux PJ, Rodseth RN. Cardiac biomarkers in the prediction of risk in the non-cardiac surgery setting. *Anaesthesia* 2014; 69: 484-493.
 3. Østervig RM, Køber L, Forberg JL, Rasmussen LS, Eugen-Olsen J, Iversen K. suPAR—a future prognostic biomarker in emergency medicine. *Scand J Trauma Resusc Emerg Med* 2015;23(Suppl 1):A31.
 4. Kjellman A, Akre O, Gustafsson O et al. Soluble urokinase plasminogen activator receptor as a prognostic marker in men participating in prostate cancer screening. *J. intern Med.* 2011 Mar;269(3):299-305.
 5. Lomholt AF, Christensen IJ, Høyer-Hansen G et al. Prognostic value of intact and cleaved forms of the urokinase plasminogen activator receptor in a retrospective study of 518 colorectal cancer patients. *Acta Oncology* 2010; 49:805-811
 6. Stephens RW1, Nielsen HJ, Christensen IJ, Thorlacius-Ussing O, Sørensen S, Danø K, Brünner N. Plasma urokinase receptor levels in patients with colorectal cancer: relationship to prognosis. *J Natl Cancer Inst.* 1999 May 19;91(10):869-74.

7. Svendsen MN, Ytting H, Brünnner N et al. Preoperative concentrations of suPAR and MBL proteins are associated with the development of pneumonia after elective surgery for colorectal cancer. *Surg Infect (Larchmt)* 2006 okt;7(5): 463-471.
8. Fidan E, Mentese A, Ozedemir F. Diagnostic and prognostic significance of CA IX and suPAR in gastric cancer. *Med oncol* 2013 jun;30(2):540
9. Fernebro E, Madsen RR, Fernö M et al. Prognostic importance of the soluble plasminogen activator receptor, suPAR, in plasma from rectal cancer patients. *Eur j cancer* 2000; 486-491.
10. Stephens RW, Nielsen HJ, Christensen IJ. Plasma urokinase receptor levels in patients with colorectal cancer: relationship to prognosis. *J natl Cancer Inst* 1999, 91, 869-874.
11. Galliera E, Drago L, Marazzi MG. Soluble urokinase-type plasminogen activator receptor (suPAR) as new biomarker of the prosthetic joint infection: correlation with inflammatory cytokines. *Clin Chim Acta* 2015 Feb 20;414:23-8. See comment in PubMed Commons below
12. Alstrup M et al. Soluble Urokinase Plasminogen Activator Receptor (suPAR) as an Added Predictor to Existing Preoperative Risk Assessments. *World J Surg.* 2019 43(3):780-790
13. Meyer J et al. suPAR is associated with risk of future acute surgery and post-operative mortality in acutely admitted medical patients. *Scand J Trauma Resusc Emerg Med.* 2018 Feb 1;26(1):11.

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